## **Coldex**<sup>®</sup>

### **INSTRUCTIONS FOR USE (EN)**

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## **Coldex**<sup>®</sup>

Instructions for use – ENGLISH

## **Contents**

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. WARNINGS
6. PRECAUTIONS
7. PRODUCT DISPOSAL
8. PRINCIPLES OF OPERATION

### **1. INTENDED USE**

Coldex<sup>®</sup> is a moist, sterile, hydrophilic wound dressing with good absorbent properties that maintains a moist wound environment in the presence of wound exudate. Coldex<sup>®</sup> has a positive and stimulating effect on the formation of healthy granulation tissue.

### **2. COMPOSITION**

Polyvinyl alcohol (PVA) foam

### **3. INDICATIONS FOR USE**

Coldex<sup>®</sup> can be used in a broad range of exuding skin defects, where primary closure is not possible or not desired including:

- Pressure ulcers
- Leg and foot ulcers
- Traumatic wounds
- Surgical wounds
- Infected wounds
- Partial thickness burns

Coldex<sup>®</sup> should only be used by medical and nursing staff.

### **4. CONTRAINDICATIONS**

Use of Coldex<sup>®</sup> is contraindicated in:

- Dry wounds
- Necrotic wounds with eschar present
- Third-degree burn wounds
- Heavily bleeding wounds

Coldex® must not be applied on epithelialized parts of a wound and the skin surrounding the wound.

## **5. WARNINGS**

- 1) Carefully monitor the patient and dressing frequently on a regular basis for signs of local or systemic infection. Make sure all wound filler pieces are removed with each dressing change to avoid the risk of retention and possible infection as a result of the residual foam.
- 2) Coldex® must not be packed too tightly or forced into any wound area.
- 3) Sometimes, after Coldex® has been applied, a slight burning sensation may appear temporarily but it will disappear spontaneously after about 15 minutes.
- 4) If after 4 weeks of treatment with Coldex® no improvement in the wound condition could be achieved overall therapy should be reassessed. First of all underlying conditions have to be detected and corrected before therapy should be continued.

## **6. PRECAUTIONS**

- 1) Do not use Coldex® if the package is breached or damaged.
- 2) Coldex® is supplied sterile; do not re-sterilize.
- 3) Do not apply alcoholic solutions on Coldex® as this compromises the integrity, composition and functionality of the foam.
- 4) Coldex® is for single use only; reuse of the device can lead to several hazards:
  - Cross contamination
  - Infection risks
  - Deterioration and alteration of the characteristics of the device due to initial use with consequent risk of injury to the patient
- 5) Coldex® is supplied moist and should be used as soon as possible after the package has been opened.
- 6) Coldex® does not contain antibiotics. It has no bacteriostatic effect on pre-existing infections, nor does it prevent the development of new infections.
- 7) Do not cut Coldex® directly over the wound cavity to avoid foam fragments from falling into the wound. Rub edges of the foam away from the open wound.
- 8) Always count the total number of dressings used to pack a wound. Document the number in the patient's records to ensure all the dressings are accounted for upon removal.
- 9) Do not use Coldex® after the expiry date.
- 10) Keep Coldex® away from sunlight and heat; store Coldex® between temperatures of 5°C and 35°C. Keep Coldex® dry during storage.

## **7. PRODUCT DISPOSAL**

In order to minimize the risk of potential infection hazards, or environmental pollution, Coldex® should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

## 8. PRINCIPLES OF OPERATION

Use clean or aseptic techniques for application according to your institutional protocol. Thorough wound assessment and cleansing should occur with each dressing change.

### Assess, cleanse and debride the wound

Assess the wound dimensions and characteristics. Debride any necrotic or eschar tissue. Cleanse the wound bed according to your institutional protocol.

Cleansing of the wound must be performed before applying the dressing for the first time and after every dressing change. Keep the skin around the wound dry.

### Dress the wound

Take Coldex® sterile out the packaging.

Cut Coldex® with sterile instruments to fit the size and shape of the wound. Do not cut Coldex® over the wound, as fragments may fall into the wound bed. Always rub Coldex® away from the wound to remove loose edges.

Gently place Coldex® in the wound and secure the material with a semi-occlusive dressing or conforming bandage e.g. by means of an elastic bandage or in the operating theatre with situation sutures at the wound edges.

Make sure Coldex® does not cover the skin around the wound or growing epithelium at the wound edge.

In deeper defects, apply several layers of Coldex®. Count number of pieces used in the wound. If multiple pieces are needed, count and record how many pieces are used to ensure all pieces will be removed with a dressing change to avoid the risk of retention and possible infection as a result of the residual foam.

Avoid over packing or forcing the dressing into the wound.

### Dressing changes

Each time the dressing is changed the wound must be thoroughly cleansed in accordance with the physician's instructions and applicable hygiene standards.

As a rule, and depending on the discharge from the wound, Coldex® has to be changed every day, but at least every two days. Coldex® should be replaced if it becomes soiled, saturated or if exudate drainage is inadequate. Where leakage occurs Coldex® should be replaced instantly.

If the wound has become clean and dry and Coldex® has turned a little rigid, Coldex® must be soaked preferably with a salt solution or sterile water before removal. Once the dressing is fully hydrated the foam can be removed.

Count number of pieces taken out of the wound. Check patient's record to see how many pieces of Coldex® were used in order to ensure all pieces are removed to avoid the risk of retention and possibility of infection caused by the residual foam.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV  
Middenweg 20  
BE-3930 Hamont-Achel

 +32 11 440 990

 1639

Email: [info@mondomed.be](mailto:info@mondomed.be)  
<http://www.mondomed.be>

	European conformity		Do not use if package is damaged
	Medical device		Consult instructions for use
	Manufacturer		Caution
	Date of manufacture		Keep away from sunlight
	Use-by date		Keep dry
	Batch code		Temperature limit
	Catalogue number		Green dot
	Sterilized using radiation		Dispose of packaging in an environmentally responsible way
	Single sterile barrier system		Composite of Paper and plastic/aluminum
	Unique device identifier		Non-corrugated fiberboard (paperboard)
	Do not re-use		Corrugated fiberboard (cardboard)
	Do not re-sterilize		

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