Fractomed[®]

Instructions for use – ENGLISH

Contents

1. INTENDED USE
2. COMPOSITION
3. INDICATIONS FOR USE
4. CONTRAINDICATIONS
5. PRECAUTIONS
6. PRODUCT DISPOSAL
7. PRINCIPLES OF OPERATION

1. INTENDED USE

Fractomed[®] is a thermoplastic splinting material used as orthoses, pre- and post-operative external immobilization devices and rehabilitation aids. Its microporosity ensures that sufficient perspiration is possible and prevents decubitus, swelling, skin inflammation, etc.

Fractomed[®] is radiolucent which means no need to remove the device to make an x-ray. Fractomed[®] is available in different sizes and forms.

2. COMPOSITION

Polyvinyl alcohol (PVA) foam

3. INDICATIONS FOR USE

Non exhaustive list

- Fractures
- Dislocations
- Sprains
- Muscle tendon contractures
- Postoperative treatment
- First aid

Fractomed[®] should only be used by medical professionals.

4. CONTRAINDICATIONS

Fractomed[®] may not be used directly on open wounds.

5. PRECAUTIONS

- 1) Do not apply alcoholic solutions on the device as this compromises the integrity, composition

and functionality of the foam.

- 2) Fractomed[®] can be remolded multiple times on the same patient.
- 3) Do not use the device after the expiry date.
- 4) Keep the device away from sunlight and heat. Keep Fractomed[®] dry during storage.

6. PRODUCT DISPOSAL

In order to minimize the risk of potential infection hazards, or environmental pollution disposal of Fractomed[®] should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards.

7. PRINCIPLES OF OPERATION

Activation of the material

Check if the steaming device is running and has reached the correct temperature (see instructions for use of the Mondomed Steaming Device).

Put the Fractomed[®] article of your choice on the grid in the steaming device and leave it in for 2 – 3 minutes. To ensure maximum modeling time, it is recommended to apply the product immediately to the patient after activation.

Make sure that the Fractomed[®] material never comes in contact with water.

Preparing the patient

Apply a (thin) stockinet on the limb to be treated. If necessary additional padding material can be used. Take the splint out of the steaming device and fit the splint to the patient. Make adjustments if necessary. The material can be cut with a plaster scissors.

After approximately 2 minutes, the splint regains its rigidity and preserves the molded shape. If needed, the splint can be put back in the steaming device to remodel it again.

Fix the splint in its position with an elastic bandage.

Any serious incident that has occurred in relation to the device should be reported to the manufacturer, who will take care for reporting the applicable Competent Authority.

 Mondomed NV
Middenweg 20
BE-3930 Hamont-Achel
Phone: +32 11 440 990
Email: info@mondomed.be

<http://www.mondomed.be>

	European conformity
	Medical device
	Manufacturer
	Use-by date
	Batch code
	Catalogue number
	Unique Device Identifier
	Single patient - multiple use
	Keep away from sunlight
	Keep dry
	Dispose of packaging in an environmentally responsible way
	Green dot

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